

# MA35

Rugged, Fully Automatic  
Infrared Moisture  
Analyzer for Use in Routine  
Operation



## Benefits

- Time-saving alternative to a drying oven – measurement results in only 10 – 20 minutes
- Simple operation
- Rugged design
- Ideal for simple routine checks
- Best price-performance ratio with Sartorius quality – Made in Germany

## Product Overview

The MA35 is the ideal basic model for quick and reliable analysis of the material moisture of liquid, pasty and solid substances using the thermogravimetric method. Sartorius's decades of experience in the area of moisture analysis guarantees quick and reliable measurement results and high user-friendliness.

# Product Information

## **Applications**

### **Ruggedness for All Areas of Application**

The compact and rugged design, combined with great precision and measurement accuracy, make the MA 35 the ideal device for uniform measuring tasks. The selection of long-lasting and rugged components, impervious to vibration and dirt, guarantees reliable measurement results, even in a harsh environment. An ideal partner for use in production, incoming goods or in the laboratory.

### **Simple Use**

Sartorius's unique, fully automatic endpoint determination makes programming a switch-off criterion unnecessary. This saves time because, on the one hand, time-consuming optimizations are no longer required, and on the other, operation of the device is simplified so that it can be used competently by anyone. This is made possible by continuous monitoring of the drying cycle by the MA35.

This monitoring then completes the measurement fully automatically using a switch-off criterion that is calculated based on the drying curve as soon as the sample weight is constant, and where despite heating, no additional loss of weight is detectable. The measurement accuracy required for this is provided by a Sartorius weighing system with readability to 1 mg, which has been specially optimized for use in high temperature ranges.

### **Quick Measurement Results**

For sample heating, the MA35 features two powerful metal tubular-shaped heating elements generating 360 watts of power. These elements facilitate quick and uniform heating of samples and thus shorten measurement times, which means a significant amount of time is saved compared to the drying oven because the finished end result can already simply be read within a few minutes.

# Technical Specifications

<b>Weighing range</b>	35 g
<b>Typical reproducibility of the weighing system</b>	From around 1 g of weight of sample taken: $\pm 0.2\%$ From around 5 g of weight of sample taken: $\pm 0.05\%$
<b>Readability</b>	1 mg; 0.01%
<b>Measurement display</b>	% moisture, % dry mass, % ATRO, g residue
<b>Typical sample quantity</b>	5 – 15 g
<b>Access to the sample compartment</b>	<ul style="list-style-type: none"><li>▪ by lifting top with wide opening angle</li><li>▪ with viewing window</li></ul>
<b>Windscreen</b>	Integrated windscreen
<b>Switch-off criterion</b>	Optional: <ul style="list-style-type: none"><li>▪ Fully automatic</li><li>▪ Manual switch-off</li><li>▪ Allotted time 0.1 – 99 min.</li></ul>
<b>Heating program</b>	Standard drying
<b>Measurement program</b>	The device is power failure-proof, so 1 program can be saved (parameter freely selectable)
<b>Temperature setting</b>	40 – 160 degrees Celsius in increments of 1 degree Celsius
<b>Sample heating</b>	Uniform with an Infrared dark radiator (metal tubular-shaped heating elements)
<b>Measured value storage</b>	Results are safely stored until the start of the next measurement
<b>Log printout</b>	<ul style="list-style-type: none"><li>▪ Non-adjustable GLP log for measurement results   Calibration   Adjustment</li><li>▪ Paper-saving brief log for measurement results</li><li>▪ An optionally available external printer is used for the printout (YDP20   YDP30)</li></ul>
<b>Inspection equipment monitoring</b>	External calibration by means of optionally available calibration weights; temperature adjustment using an optional adjustment set
<b>User guidance</b>	Through symbols that correspond to the selected functions
<b>Language selection</b>	English, German, French, Spanish, Italian, Russian
<b>Data interface</b>	RS232 C-S/V24-V28, 7 Data-bit (ASCII) 1 Parity-bit
<b>Network frequency</b>	50   60 Hz
<b>Power consumption</b>	max. 400 VA
<b>Case dimensions (W × D × H)</b>	224 × 366 × 191 mm
<b>Weight</b>	Approximately 5.8 kg

## Models Available

MA35M-230N	Standard version 230 V +15%, -20%
MA35M-115N	Standard version 115 V +15%, -20%
MA35M-115US	Standard version 115 V +15%, -20% with country-specific enhancements for the USA
MA35M-1CN230V1	Standard version 230 V +15%, -20% with country-specific enhancements for China

## Accessories

6965542	Disposable sample pans 80 units, Aluminum, round, 90 mm
6906940	Fiberglass filter for use with liquid, pasty and oily samples, 80 units
6906941	Fiberglass filter, for use with liquid and oily samples, soft Quality, 200 count, Ø 90 mm
YDS05MA	Glass panel replacement set Aluminum pane instead of glass to comply with FDA/HACCP guidelines (Conversion kit)
YDP20-OCE	Data printer
YDP30	Premium GLP laboratory printer
69Y03295	RS232 cable for connecting the printer YDP30
6906918	Color ribbon for measured value printer
6906937	Paper rolls for data printer, 5 items at 50 m
YSS43-02	External calibration weight, 30 g ± 0.3 mg DKD [Deutsche Kalibrierdienst; German Calibration Service] certificate
YTM15MA	Temperature adjustment set
YMD01B	Individual method development for the MA35 (including brief report)

### Scope of Supply


- Moisture measuring device with data interface
- Network cable
- 80 Aluminum sample pans

### Germany

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen  
Phone +49 551 308 0

### USA

Sartorius Corporation  
565 Johnson Avenue  
Bohemia, NY 11716  
Phone +1 631 254 4249  
Toll-free +1 800 635 2906

 For further Information, visit  
[www.sartorius.com](http://www.sartorius.com)